

IHE Integration Statement

Vendor	Product Name	Version	Date
TELEMIS SA	Telemis-Medical	4.2	01/05/2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow	Image Manager/Archive		
Scheduled Workflow	Image Display		
Scheduled Workflow	Performed Procedure Step Manager		
Patient Information Reconciliation	Image Manager/Archive		
Patient Information Reconciliation	Performed Procedure Step Manager		
Portable Data for Imaging	Portable Media Creator		
Portable Data for Imaging	Portable Media Importer		
Portable Data for Imaging	Print Composer		
Portable Data for Imaging	Display		
Portable Data for Imaging	Image Display		
Consistent Time	Time Client		
Mammography Image	Image Manager/Archive		
Radiation Exposure Monitoring	Image Manager/Archive		
Radiation Exposure Monitoring	Dose Information Consumer		

Internet address for vendor's IHE information : <http://www.telemis.com/ihe/>

Links to Standards Conformance Statements for the implementation

HL7	
DICOM	http://www.telemis.com/dicom

Links to general information on IHE

In North America : http://www.ihe.net	In Europe : http://www.ihe-europe.org	In Japan : http://www.jira-net.or.jp/ihe-j
--	---	--