

IHE Integration Statement

Vendor	Product Name	Version	Date
Telemis SA	TM-Medical	4.1	May 12, 2010

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Manager/Archive	
Consistent Presentation of Images	Print Composer	
Consistent Time	Time Client	
Evidence Documents	Image Display	
	Image Manager/Archive	
Key Image Note	Image Display	
	Image Manager/Archive	
Mammography Image	Image Manager/Archive	
Patient Information Reconciliation	Image Manager/Archive	
	Performed Procedure Step Manager	
Portable Data for Imaging	Display	
	Image Display	
	Portable Media Creator	
	Portable Media Importer	
	Print Composer	
	Report Reader	
Radiation Exposure Monitoring	Dose Information Consumer	
	Image Manager/Archive	
Scheduled Workflow	Image Display	
	Image Manager/Archive	
	Performed Procedure Step Manager	
Simple Image and Numeric Report	Report Reader	
	Report Repository	

Internet address for vendor's IHE information: <http://www.telemis.com/content/default.asp?page=6251>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	http://www.telemis.com/dicom

Links to general information on IHE

In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org
--	---	--