

IHE Integration Statement		Date: 01/05/2008
Vendor	Product	Name Version
TELEMIS SA	TELEMIS MEDICAL	4.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	Image Manager/Archive	None
Scheduled Workflow	Image Display	None
Scheduled Workflow	Performed Procedure Step Manager	None
Patient Information Reconciliation	Image Manager/Archive	None
Patient Information Reconciliation	Performed Procedure Step Manager	None
Access to Radiology Information	Image Manager/Archive	None
Key Image Note	Image Manager/Archive	None
Key Image Note	Image Display	None
Evidence Documents	Image Manager/Archive	None
Evidence Documents	Image Display	None
Portable Data for Imaging	Image Display	None
Portable Data for Imaging	Portable Media Creator	None
Portable Data for Imaging	Portable Media Importer	None
Portable Data for Imaging	Print Composer	None
Portable Data for Imaging	Report Reader	None
Portable Data for Imaging	Display	None
Consistent Presentation of Images	Print Composer	None
Simple Image and Numeric Report	Report Reader	None
Simple Image and Numeric Report	Report Repository	None
Consistent Time	Time Client	None
Mammography Image	Image Manager/Archive	None
Internet address for vendor's IHE information : http://www.telemis.com/ihe/		
Links to Standards Conformance Statements for the implementation		
DICOM:	http://www.telemis.com/dicom	
Links to general information on IHE		
In North America : http://www.ihe.net	In Europe : http://www.ihe-europe.org	In Japan : http://www.jira-net.or.jp/ihe-j