



# **Telemis-Medical**

## IHE Integration Statement

Version 4.8

Revised June 2017

**Telemis s.a.**  
Avenue Athena, 2  
B-1348 Louvain-la-Neuve  
Belgium

**Tel.** +32-10-48 00 11  
**Mail** [info@telemis.com](mailto:info@telemis.com)  
**Web site** <http://www.telemis.com>



Extending Human Life...



## About this document

Copyright © 1999 - 2017 Telemis s.a., All rights reserved.

### Disclaimer

Telemis believes that the information in this document is accurate at the release date but shall not be liable for errors contained herein or for any loss or damages in connection with this document.

Telemis reserves the right to revise and to make changes to this document at any time, without any obligation to notify of such revisions and changes. The latest up-to-date version is always available at [www.telemis.com/ihe](http://www.telemis.com/ihe).

### Contact information

Telemis s.a.  
Avenue Athena, 2  
B-1348 Louvain-la-Neuve  
Belgium

**Tel.** +32-10-48 00 11

**Mail** [info@telemis.com](mailto:info@telemis.com)

**Web site** <http://www.telemis.com>



## TABLE OF CONTENTS

ABOUT THIS DOCUMENT.....	2
1 IHE INTEGRATION STATEMENT .....	4



# 1 IHE Integration Statement

Vendor	Product Name	Version	Date
Telemis s.a.	TM-Medical	4.8	01/06/2017

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration profiles	Actors	Options	Last Connectathon Tested
Consistent time	Time client	None	EU 2017
Key Image Note	Image Display	None	EU 2017
Mammography Image	Image Manager/Archive	None	EU 2017
Access to Radiology Information	Image Manager/archive	None	EU 2009
Portable Data for Imaging	Image Display	None	EU 2016
	Portable Media Creator	None	EU 2017
	Portable Media Importer	None	EU 2017
	Display	None	EU 2011
	Print composer	None	EU 2011
CT/MR Perfusion with Imaging Contrast Image	Manager/Archive	None	EU 2016
Patient Information Reconciliation	Image Manager/Archive	None	EU 2013
	Performed Procedure Step Manager	None	EU 2013
Radiation Exposure Monitoring	Dose Information Consumer	None	EU 2017
	Image Manager/Archive	None	EU 2017
Scheduled Workflow	Image Display	None	EU 2013
	Image Manager/Archive	None	EU 2013
	Performed Procedure Step Manager	None	EU 2013
Scheduled Workflow.b	Image display	None	EU 2017
	Image Manager/Archive	None	EU 2017
	Performed	None	EU 2017



	Procedure Step Manager		
Cardiology Evidence Documents	Image display	None	EU 2016
	Image display	Echo Evidence	EU 2016
Echocardiography Workflow	Image display	None	EU 2016

Internet address for vendors IHE information: <a href="http://www.telemis.com/">http://www.telemis.com/</a>	
Links to Standards Conformance Statements for the Implementation	
<b>HL7</b>	
<b>DICOM</b>	<a href="http://www.telemis.com/ihe">http://www.telemis.com/ihe</a>
Link to general information on IHE: <a href="http://www.ihe.net">http://www.ihe.net</a>	

